VOCO

Ufi Gel hard C

GB

Instructions for use

CE**Product description:**

Ufi Gel hard C is a cold-curing permanently hard relining material for denture relinings. It is simple and quick to use for direct as well as indirect relinings. **Ufi Gel hard C** contains no methyl methacrylates, thus minimising the risk of allergies and irritation of the mucosa. The set additionally contains an adhesive to achieve a permanent and stable bond.

Indications:

- hard, permanent total or partial relinings to restore the functions of partial and complete dentures
- lengthening of denture margins

Application:**1. Preparation of the denture**

Check occlusion and carry out corrections if necessary. Clean denture thoroughly with a brush and dry. Roughen all areas to be relined incl. buccal and labial margins with a suitable bur. Then clean and dry. Insulate all areas not to be relined (e.g. artificial teeth, labial and buccal areas beyond the roughened surfaces), e.g. with vaseline, to prevent unwanted contact with **Ufi Gel hard C**. **Partial clasp dentures, anchors and attachments of dentures have to be blocked out with silicone (liquid) or wax with a low melting point.**

2. Application of the adhesive

Exchange the transportation cap against the brush cap. Apply the adhesive to all surfaces to be relined and let it air-dry (30 s). **Attention: Close the adhesive bottle tightly immediately after use and store in an upright position.**

Important: For a secure bond use only Ufi Gel hard C adhesive.

3. Preparation of the cartridge

Put the **Ufi Gel hard C** cartridge into a **VOCO** dispenser (type 2) or any extruder gun with a suitable plunger. Slight differences in the filling levels of the two cartridge chambers may occur for technical reasons.

Therefore, remove the transportation cap from the cartridge and extrude material until both chambers deliver uniformly. Then fix and lock the mixing tip with a 90° turn. The material is extruded and automatically mixed in the correct ratio by pressing the dispenser lever. After use leave the mixing tip on the cartridge or refit the transportation cap.

4. Application of Ufi Gel hard C

Apply **Ufi Gel hard C** evenly onto the prepared margin and/or base of the denture, avoid leaving excess material and remove any with a suitable instrument. Re-insert denture and have the patient exert slight pressure in occlusal position for 1 min and then carry out functional, chewing and swallowing movements for 2 min. Make sure that no material flows into the throat when relining an upper denture, especially at the transition from hard to soft palate (A-line).

4.1 Partial or total dentures with undercuts

Remove denture 3.5 min after the beginning of application and remove excess material immediately with scissors or a scalpel. Re-insert into the mouth for a final occlusion, check and allow to cure for 2 min approx. Curing can also be completed in warm water, e.g. in a pressure pot at approx. 40 °C, if the relined denture is not to be re-inserted. **Do not let material cure in contact with air** since oxygen will inhibit complete curing of the surface.

4.2 Total dentures without undercuts

Excess material can be removed intra- or extraorally before final polymerisation. Intraorally: remove excess material after 3 min (from beginning of application) with a suitable instrument. Let the denture cure for a further 2 min in the mouth until **Ufi Gel hard C** is completely cured. Extraorally: see position 4.1

5. Finishing

The whole of the relined surface (also basal) is to be treated using standard tools for finishing (tungsten carbide bur, silicone polisher) and polishing (polishing paste, polishing disc).

beginning of application	0,5 min	apply
insert denture into the mouth	1 min	exert pressure
	2 min	functional movements
remove denture if necessary	1,5 min	remove excess
insert denture into the mouth	2 min	curing
remove denture	2,5 min	finishing, polishing
	9,5 min	

extraorally intraorally intra- or extraorally

Information, precautionary measures:

- Store adhesive bottle tightly closed and in an upright position.
- Avoid contact of adhesive with skin. Wash contact points thoroughly with water and soap.
- In case of contact with the eyes rinse thoroughly with water and consult an ophthalmologist.
- **Ufi Gel hard C** contains benzoyl peroxide (catalyst), hydroxyethylmethacrylate and acetone (adhesive).
- **Ufi Gel hard C** should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- The stated time periods refer to a room temperature of 25 °C and a temperature in the mouth of 35 °C. Polymerisation will be slower at lower temperatures, higher temperatures will accelerate polymerisation.
- Do not sterilise mixing tips in ultrasonic bath.

Storage:

Store at temperatures of 4 °C to 23 °C. Do not use after expiry date.

Cleaning indications:

Dentures relined with **Ufi Gel hard C** can be cleaned with the usual cleansers and procedures. This applies to domestic as well as professional cleaning.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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VC 60 BB2215 E1 0716 99 © by VOCO
VOCO GmbH