

# **Identium®** Light



Here used in combination with Identium® Heavy in the double mix impression technique

## **Identium**® – the best of two worlds.

Identium® Light is a low viscosity, addition-curing elastomeric impression material based on Vinylsiloxanether®.

Identium<sup>®</sup> Light has been specially developed for the one-step impression technique and combines the best of two worlds: all the advantages of polyether materials coupled with those of A-silicones.

The highly elastic properties of Identium<sup>®</sup> Light ensure dimensionally accurate recovery and very easy removal. Identium<sup>®</sup> Light is also odorless and has a neutral taste.

Based on the entirely new material of Vinylsiloxanether® Identium® Light provides the following benefits:

• Excellent flowability combined with exceptional hydrophilicity (lowest achievable contact angle of less than 10° after 1 second), ensuring optimal flow behavior in a moist environment, even in the narrowest sulcus crevices.

### Vinylsiloxanether® (VSXE) Precision impression material

### SUGGESTED IMPRESSION TECHNIQUE

One-step impression technique (simultaneous technique)\*

### **INDICATIONS**

Implant impressions

Impressions for crowns/bridges/inlays/onlays and veneer preparations

Impressions about full or partial denture

highly recommended recommended

\*in combination with Identium® Heavy




### **Precision impression material (VinyIsiloxanether®)**

### The Identium® Light Advantages and Benefits

#### **Product Advantages**

- Excellent flowability
- Exceptional hydrophilicity
- Optimal elastomeric properties
- Balanced setting behavior
- Odorless and neutral taste
- Accurate model fabrication
- Easy handling

# Contact angle measurement




Identium® after 1 sec.: contact angle < 10°

Identium® Light is a low viscosity impression material available in the 1:1 cartridge. It is particularly suitable for the double mix impression technique (in combination with high viscosity Identium® Heavy), for taking impressions for the fabrication of crowns/bridges, inlays/onlays and veneers, for reline impressions as well as impressions for removable restorations.

#### **Benefits to the Practitioner**

- Material flows into the narrowest sulcus crevices with high positional stability, constantly guaranteed throughout the entire working time.
- Optimal flow behavior in a moist environment, lowest achievable contact angle (less than 10° after 1 second!).
- Dimensionally accurate recovery and easy removal from the patient's mouth.
- Long working time with short intraoral setting time (double-snap effect).
- More comfortable for the patient: no gag reflex, no erratic movement.
- Minimum risk of fracture when removing the impression.
- Compatible with commercially available disinfectants (no swelling), can be poured immediately.

### Identium<sup>®</sup>. A perfectly balanced impression material.



### **Technical Data**

Product name	Mixing ratio	Working time at 23° C/ 74° F ≤	Working time at 35° C/ 95° F ≤	Intraoral setting time at 35° C/ 95° F ≥	Total setting time* ≥	Hardness (approx.) Shore	Linear dimensio- nal change (max.) %	Elastic Recovery test (approx.)	Strain in compres- sion (approx.)
ldentium® Light	1:1	2 minutes	80 seconds	3 minutes 30 seconds	5 minutes 30 seconds	A 46	- 0.20	99.0	3.8

Product name	ISO 4823	Consistency	Mixing ratio and total	Mixing technique	Mixing element	Color	Odor/taste
		(approx.)	content volume per				
		mm	unit				
Identium® Light	Type 3,	42	1:1, 50 ml cartridge	Applyfix® 4 dispensing	Yellow mixing tip	Durnlo	Neutral
iuciiiiuiii Ligiii	light-bodied			gun DS-50 1:1/2:1	MB Ø 4.2 mm	i uipie	

<sup>\*</sup> Total setting time (removal from the mouth) from the start of mixing.