


MAP One

System for
the professional
placement
of endodontic
repair materials



New!
Simplified
Economic

MAP One is an ultra simplified version of the original MAP System with only one syringe and one «Memory Shape» needle for all your clinical repair cases.



Swiss quality
dental products

MAP One

The growing demand from worldwide general practitioners for the MAP System encouraged us to offer a new version of the product, more affordable and better adapted to the needs of the greatest number. Therefore, MAP One is a simplified starting kit bringing together the required articles enabling any dentist to place endodontic repair materials with precision and no waste. As a reminder, the MAP System provides a unique and efficient method for placing such materials (e.g. MTA, cements, bioactive dentine substitutes, calcium hydroxide...) by orthograde or retrograde obturation.



MAP One Benefits

- Optimal usage of the material (no waste)
- High efficiency and accuracy placement
- Compatible syringe with all MAP System needles
- Memory shape NiTi needle bendable to any required curvature and clinical cases

MAP One

REF
20296

1 stainless steel syringe
1 NiTi «Memory Shape» needle N°1 (rouge / ext. Ø1.10 mm)
8 plastic plungers
1 cleaning curette

Refill plastic plungers (16 pieces)

20201



Produits Dentaires SA
Vevey . Switzerland

MAP SYSTEM
MICRO - APICAL PLACEMENT

www.pdsa.ch
© PD 160809