

Temporary Fillingmaterial



Economical & reliable



coltosol® F

Simple Application, reliable effect

With Fluoride No Eugenol

Easy application

- Creamy consistency for easy adaptation in the cavity
- Good retention in the cavity
- No mixing, no waste of material
- Simple setting by the action of saliva
- Easy removal in large portions
- Long shelf life



- Fluoride release
- Radioopacity for the perfect control of the treatment result
- Perfect marginal seal
- Eugenol-free, no effect on the curing of composites
- Ideal surface hardness for provisional fillings in posteriors as well as for root canal sealings

The choice is yours

The required quantity of Coltosol F is either dosed with the syringe or taken directly from the glass jar with an instrument. Quick, easy and clean.



www.coltene.com © 2014 Coltène/Whaledent AG





Instant Application

After rolling the dispensed Coltosol F into a ball with your fingers, simply press it into the prepared cavity.



Easy removal

Coltosol F can be easily removed in large pieces using a bur or an excavator without leaving any residues.

Economical: Optimal value for money



Order Information:

Coltosol F in syringes:

5930 Coltosol F

5 syringes of 8 g each

Coltosol F in glass jar:

5906 Coltosol F Single pack

1 glass jar of 38 g

5911 Coltosol F Eco Pack

3 jars of 38 g each

Coltène/Whaledent AG

Feldwiesenstrasse 20 9450 Altstätten / Switzerland Tel +41 71 757 5300 Fax +41 71 757 5301 info.ch@coltene.com

Coltène/Whaledent Inc.

235 Ascot Parkway Cuyahoga Falls, Ohio 44223 / USA Tel +1 330 916 8800 Fax +1 330 916 7077 info.us@coltene.com

Coltène/Whaledent Ltd.

The President Suite Kendal House, Victoria Way Burgess Hill, West Sussex RH15 9NF / UK Tel +44 1444 235486 Fax +44 1444 870640 info.uk@coltene.com

